



Indirect Bonding Tray

Used for printing orthodontic transfer trays for precise positioning during bracket placement. Drastically reduce chair time and increase patient comfort during the orthodontic bracket setting process.

Available Materials :

KeyOrtho IBT® by Keystone Industries is a biocompatible, tasteless, and odorless material that combines the strength and precision to accurately set brackets, with the flexibility and non-stick formula needed for easy release.

System	Origin® One Dental
Material	KeyOrtho IBT® by Keystone Industries
Build Time	25m
Material Used	46.31ml
Built plate include	22 models
Print Up to	465 Parts/Day*

*10h/day+ 1 overnight tray print, 95% utilization, 261 workdays/year

ISO 9001:2015 Certified

© 2022 Stratasys. All rights reserved. Stratasys, the Stratasys Signet logo and Origin one Dental are registered trademarks of Stratasys Inc. All other trademarks are the property of their respective owners, and Stratasys assumes no responsibility with regard to the selection, performance, or use of these non-Stratasys products. Product specifications are subject to change without notice

STRATASYS.COM