



Soft Splints

The use of occlusal appliances is quickly becoming the standard of care. Flexible splints, night guards, and bleaching trays printed with the Origin One Dental provide a comfortable fit, flexibility to engage undercuts, and excellent resistance to wear.

Available Materials :

KeySplint Soft Clear® by Keystone Industries is strong, flexible, fully biocompatible (510K Cleared Class IIa Certified), and easy to polish and clean.

System	Origin® One Dental
Material	KeySplint Soft Clear® by Keystone Industries
Build Time	1h 58m
Material Used	39.6ml
Built plate include	9 models
Print Up to	51 Parts/Day*

*10h/day+ 1 overnight tray print, 95% utilization, 261 workdays/year

ISO 9001:2015 Certified

© 2022 Stratasys. All rights reserved. Stratasys, the Stratasys Signet logo and Origin one Dental are registered trademarks of Stratasys Inc. All other trademarks are the property of their respective owners, and Stratasys assumes no responsibility with regard to the selection, performance, or use of these non-Stratasys products. Product specifications are subject to change without notice.

STRATASYS.COM

