



Hard Splints

Fabricate rigid dental splints for bruxism therapy, retainers, and other cases where tooth immobilization is needed.

Available Materials :

KeySplint Hard® By Keystone Industries is a fully biocompatible, easy-to-polish and clean, strong, and abrasion-resistant resin.

System	Origin® One Dental
Material	KeySplint Hard® by Keystone Industries
Build Time	1h 45m
Material Used	49.74ml
Built plate include	9 models
Print Up to	56 Parts/Day*

*10h/day+ 1 overnight tray print, 95% utilization, 261 workdays/year

ISO 9001:2015 Certified

© 2022 Stratasys. All rights reserved. Stratasys, the Stratasys Signet logo and Origin one Dental are registered trademarks of Stratasys Inc. All other trademarks are the property of their respective owners, and Stratasys assumes no responsibility with regard to the selection, performance, or use of these non-Stratasys products. Product specifications are subject to change without notice.

STRATASYS.COM

