

KeyGuide®

KeyGuide is ideal for fabricating transparent and stable surgical guides, which enable the placement of implants at a precise angle and depth.



BIOCOMPATIBLE



STRONG



EASY TO POLISH



AUTOCALVABLE


Class I

#4200001 (0.5kg bottle)

#4220001 (1kg bottle)




FLEXURAL MODULUS (MPa)



✓ AVAILABLE IN THE USA, CANADA, EU & AUS

FLEXURAL STRENGTH (MPa)



keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
 EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
 EU: +49 77 31 91 21 01