



# TrueDent™ Denture

Monolithic, polychromatic digital denture printed using TrueDent resins.

TrueDent is an FDA cleared (Class II) resin developed for 3D printing of dentures and temporary crowns and bridges on the J5 DentaJet platform. It enables batch production of highly aesthetic, monolithic, multicolor parts on a single mixed-part, high-capacity tray.

# TrueDent™ Denture

<b>System</b>	Stratasys J5 DentaJet
<b>Build Time</b>	High Quality Speed mode: 2 hours, 43 mins
<b>Model Materials Used</b>	TrueDent™ Cyan, TrueDent™ Magenta, TrueDent™ Yellow, TrueDent™ White, TrueDent™ Clear
<b>Support Material Used</b>	TrueDent™ Support

FDA cleared (Class II)  
ISO 20795-1:2013 Certified  
ISO 10477:2020 Certified

© 2023 Stratasys. All rights reserved. Stratasys, Stratasys signet, J5 DentaJet, TrueDent Cyan, TrueDent Magenta, TrueDent Yellow, TrueDent White, TrueDent Clear are trademarks of Stratasys Ltd. and/or its subsidiaries or affiliates and may be registered in certain jurisdictions. All other trademarks belong to their respective owners.  
PC\_PJ\_Truedent Denture\_A5\_0123a

STRATASYS.COM

